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## **1. Introduction**

The purpose of this document is to provide detailed instructions for health care facility reprocessing of the intramedullary reamers marketed by MPS Precimed.

It also contains instructions for cleaning, inspection, sterilisation, storage and use of the instruments.

The cleaning and sterilisation processes outlined in this leaflet have been validated by MPS Precimed. Other reprocessing methods may be effective, but they must be validated by the end-user under its responsibility before being used. Furthermore, the end-user must comply with the laws and regulations of countries with stricter reprocessing requirements than those specified in this leaflet.

---

## 2. Warnings

A complete reprocessing is necessary prior to each use, and this also applies to new devices.

The instruments are not intended for use in the treatment of children or expectant mothers. Nevertheless, within the strict framework surgery of limb lengthening with no aesthetics purpose, it is allowed to use the instruments on children over the age of 8 years and where the growth plates are joined.

The use of personal protective equipment is strongly recommended when handling contaminated or potentially contaminated instruments.

The maximum permissible exposure temperature for the instruments is 137°C.

Highly alkaline ( $\text{pH} > 11$ ) and hypochlorous solutions must be avoided as they promote corrosion of metallic parts.

The use of metal brushes during cleaning is strictly forbidden as this will damage the surface of the instruments.

Complex instruments with long, narrow tubes and blind holes require special attention during cleaning.

Instruments with dried-on soil are more difficult to clean.

The instruments are designed for a total number of 75 uses including the necessary reprocessing (cleaning, disinfection, sterilisation). When an instrument has reached the end of its useful life, it is advisable to clean it of all biological substances and to discard it in accordance with the prevailing laws and regulations.

Instruments which show signs of wear or damage must not be used. In particular, if cracks or scratches are visible on the titanium-nickel alloy (Nitinol) shaft, the instrument must not be used under any circumstances and must be absolutely replaced by a new one. Tests prove that these defects are the cause of shaft breakage, which may take place during use in the patient's medullary cavity with the risk of not being able to remove all the fragments.

A guide wire is necessary to ensure the instruments are properly used.

---

The reamer head must be completely introduced in the canal before placed in action to avoid an excessive reaming of the anterior cortex.

An excessive torque and/or an abrupt stoppage of the instrument can cause a fracture of the Nitinol shaft and its fragmentation.

The coupling system between the reamer head and the flexible drive shaft for the modular version is specific to MPS Precimed and is not compatible with other modular systems. Therefore, it is strictly prohibited to use a competitor's flexible drive shaft with an MPS Precimed reaming head, and vice versa.

### **3. Reprocessing at the point of use**

Remove excess soil with lint-free single-use wipes as soon as possible after use.

Transport the instruments in moist paper and avoid jolts and mechanical damage.

## 4. Cleaning

The instruments must be cleaned as soon as possible using one of the two techniques specified below in order to prevent the soil from drying.

Only agents, solutions and detergents proven to be effective may be used. The agents, solutions and detergents specified below are those used by MPS Precimed in validating these instructions.

### a. Manual cleaning

#### Equipment and consumables:

- Alkazyme® enzymatic solution
- Neodisher® MediClean alkaline detergent
- Osmosed water or the chemical and microbiological equivalent
- Soft nylon bristle brush
- Single-use towelettes
- Ultrasonic bath

---

### Method:

- Rinse and brush the instrument for at least 90 seconds (see details in the appendix on page 93) to remove excess soil while completely submerging it in a bath of room temperature (15 to 25°C) osmosed water or the chemical and microbiological equivalent.
- Completely submerge the instrument in the Alkazyme® enzymatic solution at 1% volume per volume (prepared according to the manufacturer's instructions) for 18 minutes at room temperature (15 to 25°C, ideally 20°C). Brush all the surfaces of the instrument (see details in the appendix on page 93) with a soft nylon bristle brush for at least 90 seconds until visible soil is removed. Make sure the tubes and holes are thoroughly cleaned.
- Remove the instrument from the solution and rinse in room temperature (15 to 25°C) running water at for at least 3 minutes. Submerge the instrument in an ultrasonic bath of MediClean detergent at 1% volume per volume (prepared according to the manufacturer's instructions) for 18 minutes at a maximum temperature of 40°C.

- Rinse the instrument with room temperature (15 to 25°C) osmosed water for at least 3 minutes by first pouring into the cannulation 3 volumes of 50 mL running water each. Ensure that the water enters the tubes and that the blind holes are filled and emptied several times.
- Thoroughly dry the instrument with single-use towelettes. Visually inspect the instrument in a well-lit area to verify that no soil remains on the device.
- Repeat the previous 5 steps if there is persistent visible soil.

#### b. Automated cleaning

##### Equipment and consumables:

- Alkazyne<sup>®</sup> enzymatic solution
- Neodisher<sup>®</sup> MediClean alkaline detergent
- Osmosed water and purified water or the chemical and microbiological equivalent of each
- Soft nylon bristle brush
- Ultrasonic bath
- Validated washer-disinfector serviced in accordance with prevailing local procedures

---

Method:

- Rinse and brush the instrument for at least 30 seconds to remove excess soil while completely submerging it in a bath of room temperature (15 to 25°C) osmosed water or the chemical and microbiological equivalent.
- Submerge the instrument in an ultrasonic bath of Alkazyme® enzymatic solution at 0.5% volume per volume (prepared according to the manufacturer's instructions) for 15 minutes at room temperature (15 to 25°C). Brush all surfaces of the instrument with a soft nylon bristle brush for at least 30 seconds in the solution until visible soil is removed. Make sure the tubes and holes are thoroughly cleaned.
- Thoroughly rinse the instrument with purified water at room temperature (15 to 25°C) for at least 1 minute.
- Load the instrument into the washer-disinfector, and arranging it so that tubes and holes can drain.
- Run a washer-disinfector cycle for 10 minutes at temperature of 93°C with MediClean detergent at 0.5% volume per volume (prepared according to the manufacturer's instructions)
- When unloading, visually inspect the tubes, holes and all other areas of the instrument which are difficult to access in a well-lit area to verify that all visible soil has been removed. If necessary, repeat the cycle and/or clean manually.

- Make sure that the instrument is completely dry. If necessary, use single-use towelettes to remove any remaining traces of moisture.

## **5. Inspection**

Visually inspect the instruments in a well-lit area to detect any signs of corrosion, damage or wear. The cutting edges must be free of nicks and present a continuous edge.

Remove the damaged instruments or instruments with dull cutting edges, clean them of all biological substances and discard them in accordance with the prevailing laws and regulations.

## **6. Sterilisation**

### **a. Preparation**

The instruments are placed in sterilization trays or in specialized instrument trays and then wrapped in medical quality packaging for standard steam sterilization, according to the double wrapping technique. Also be sure as far as possible not to allow the instruments to collide against each other.

## b. Cycles

The instruments must be sterilised by moist steam autoclave using a pre-vacuum cycle (ISO 17665-1). The autoclave must be validated, serviced and calibrated in accordance with local prevailing procedures.

The following cycles have been validated to ensure a sterility assurance level (SAL) of  $10^{-6}$ :

<b>Cycle type</b>	<b>Temperature</b>	<b>Steam exposure time (minimum)</b>	<b>Drying time (minimum)</b>
Pre-vacuum	132°C	4 minutes	20 minutes
	134°C	3 minutes	20 minutes
	134°C	18 minutes *	20 minutes

\*Steam sterilisation parameters recommended by the World Health Organisation (WHO) for instruments subject to risk of TSE/CJD (transmissible spongiform encephalopathy and Creutzfeldt-Jakob disease) contamination.

## **7. Storage**

Store the reprocessed instruments away from dust, mildew, insects and vermin, and away from extremes of temperature and humidity.

The maximum allowable storage time between uses for sterile instruments must be defined by the health care facility.

## **8. Use**

MPS Precimed instruments are to be used only by qualified staff fully trained to use them in surgery and for appropriate surgical procedures.

The instruments are intended to be used on adult patients.

It is necessary to put a guide wire of a suitable diameter in the medullary cavity, which allows guiding the reamer and locking the reamer head on the flexible drive shaft for the modular version. The guide wire's diameter is 2 mm for the humerus and 3 mm for the femur or the tibia.

---

During the surgical procedure, check frequently the guide wire's position using a image intensifier to prevent the guide wire from advancing and penetrating unintentionally into the neighbouring tissues. Check also the advance of the reamer head with the image intensifier.

The medullary cavity is initially reamed with the smallest diameter and then the diameter is increased in 0.5 mm steps up to the desired size.

During the use, bone debris can accumulate in the cannulation ; to avoid this, it is recommended to clean the cannulation.

In case of resistance, do not force the instrument, but reverse the drive direction and remove it from the medullary cavity.

The surgical power tool will be connected to the proximal end of the intramedullary reamer by a coupling (Big AO, Stryker / Zimmer-Hall, Hudson).

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