FRAISES – Notice d'utilisation
ФРЕЗЕР – Инструкции за употреба
FRÉZY – Návod k použití
FRÆSERE – Brugsanvisning
FRÄSEN – Gebrauchsanleitung
ΦΡΑΟΥΛΕΣ – Οδηγίες χρήσης
REAMERS – Instructions for use
FRESAS – Manual de uso
RIIMERID – Kasutusjuhend
LYUKTÁGÍTÓK – Használati útmutató
FRESE – Istruzioni per l'uso
FREZOS – Naudojimo instrukcijos
FRĒZES – Lietošanas rokasgrāmata
BOR – Bruksanvisning
WIERTŁA – Instrukcja użytkowania
REAMERS – Manual de Utilização
FREZE – Instrucțiuni de utilizare
VÝSTRUŽNÍKY – Návod na použitie
POVRTALA – Navodila za uporabo
FRÄS – Bruksanvisning
FREZELER – Kullanım Kılavuzu

FR
BG
CS
DA
DE
EL
EN
ES
ET
HU
IT
LT
LV
NO
PL
PT
RO
SK
SL
SV
TR

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REAMERS – Instructions for use

1. INTRODUCTION	57
2. OVERVIEW	57
3. WARNINGS AND PRECAUTIONS FOR USE	58
4. INSTRUCTIONS FOR USE	59
5. REPROCESSING	59
6. STORAGE AND DISPOSAL	62
7. DEFINITION OF SYMBOLS	62
8. SERIOUS INCIDENTS	63

EN

1. INTRODUCTION

This document provides instructions for use and reprocessing for users of the acetabular reamers marketed under the name of the manufacturer.

It contains the instructions for use of the instruments and also instructions for cleaning, inspection, sterilisation and storage.

The cleaning and sterilisation procedures described in this manual have been validated by the manufacturer. Other reprocessing methods may be suitable, but they must be validated beforehand by the end user and remain the responsibility of the latter. Moreover the end user must comply with the laws and regulations in countries with stricter reprocessing requirements than those specified in this manual.

2. OVERVIEW

2.1 Description and composition

Acetabular reamers are reusable instruments designed to prepare the site for placement of an acetabular cup during orthopaedic surgery.

The acetabular reamers are composed of stainless steel.

The principal action is obtained by cutting edges arranged around the rim of the hemispheric device, their shape enabling debris to be collected inside the latter.

2.2 Intended use

The instruments described in this manual are used to prepare the acetabulum to fit the shape of the acetabular implant via a reamer handle connected to a surgical motor providing rotational power.

They are designed to reshape the hip bone, creating a cavity into which the acetabular implant

will be fitted during orthopaedic surgery in adult patients.

The instruments have been validated for 60 cycles of use, including the required reprocessing (cleaning and sterilisation), without encountering any functional problem. Any further use is the responsibility of the user. The use of the instrument is not limited by the number of validation cycles, but must only be reused if inspection after reprocessing shows the absence of signs of wear or damage such as described in chapter 5.3.

They must only be used by qualified personnel fully trained in their use and in the appropriate surgical procedures in an operating theatre.

2.3 Contra-indications

Acetabular reamers are not intended for use on children because reaming the growth plate impairs bone growth. The use of acetabular reamers in pregnant women is also forbidden because of the contra-indication to anaesthesia in this patient population.

Acetabular reamers are not intended for use in maintenance operations such as for sharpening cutting edges or straightening teeth.

2.4 Performance characteristics

The cutting diameter of acetabular reamers is from 36 to 80 mm. They will withstand a maximum axial load of 250N and 17Nm torque.

The instrument is designed for axial load only, heavy lateral load should not be applied on the reamer handle during reaming as this could lead to premature wear or even breakage.

3. WARNINGS AND PRECAUTIONS

FOR USE

3.1 Prerequisites before use

A full cycle of reprocessing consisting of cleaning, inspection and sterilisation must be performed before each use and this applies to new instruments as well. The use of non-reprocessed instruments can damage the patient's health by causing infection.

An inspection must be performed before each use. Any instrument showing signs of wear or damage must under no circumstances be used and must be replaced by a new instrument.

Before use, the acetabular reamers, the connections between the various parts, including that with the surgical motor, must be duly checked. If the coupling between the parts breaks during surgery, the duration of the procedure will be extended.

3.2 Precautions for use

The use of personal protective equipment is strongly recommended during the handling of contaminated or potentially contaminated instruments.

Careful handling of the acetabular reamers is recommended to avoid the risk of injury and tearing of surgical gloves on the cutting edges.

Avoid contact between instruments and any other equipment or tools present, specifically within the location around the open wound during use.

Attention should be paid to the retractors which are particularly exposed during an anterior minimally invasive approach as this could lead to the creation of metal particles entering the surgical cavity due to friction and shock.

The use of a hammer or any other tool to assist

the reaming process should be avoided at all times as this can cause severe damage to the instrumentation.

Care must be taken when removing the acetabular reamer from the surgical cavity, as the retractors may obstruct the reamer passage and lead to premature wear or even breakage of the devices.

3.3 Environment for use

The products covered by this manual are intended for use in operating theatres in sterile areas, in the human body including contact with vital fluids such as blood.

During the entire life cycle of the product, the maximum permissible temperature for instruments is 137°C, beyond which they may deteriorate. In addition, strongly alkaline solutions (pH > 11) and hypochlorite solutions should be avoided because they promote corrosion of metallic parts.

3.4 Reprocessing precautions

During cleaning processes, the use of personal protective equipment is strongly recommended because contact with instruments soiled with human blood may cause infection and contamination of staff.

It is strictly prohibited to use metal brushes for cleaning, they cause premature wear of the instrument.

Acetabular reamers are simple instruments which require no special attention during cleaning.

The cleaning process of the acetabular reamers must be started as soon as possible after use, soiled instruments which have dried are more difficult to clean.

As far as possible, any abrupt contact with the

cutting parts of the acetabular reamers must be avoided because there is a risk of injury and tearing of surgical gloves.

4. INSTRUCTIONS FOR USE

Before use, consult the precautions for use described in chapter 3.2.

The good condition of the acetabular reamer must be checked before use.

All the instruments must be used and connected to a reamer handle with a crossbar connexion, in a sterile area of the operating theatre.

The reaming procedure on the acetabulum shall be performed in compliance with current surgical techniques.

Before reaming ensure the acetabulum is accessible and allowed sufficient femoral movement. Particular attention should be paid to ensure that the femoral capsule is released and the femoral head removed. In this configuration the femur could be easily moved providing better access.

It is recommended to start the power tool gradually before applying axial force. This will reduce the risk of blocking the reamer in the bone due to excessive torque, even if the bone composition of the acetabulum is dense. This also allows early detection and reaction in case of interference.

During reaming with the acetabular reamer, it is recommended to use a sweeping front to back movement to obtain a smooth surface.

5. REPROCESSING

Before reprocessing, consult the reprocessing precautions described in chapter 3.4.

5.1 Reprocessing at the point of use

It is recommended to remove excess soil with lint-free disposable wipes as promptly as possible after the end of the procedure.

Transport the instruments wrapped in damp paper between the various sites of use and/or reprocessing, avoiding shocks and mechanical damage.

5.2 Cleaning

The instruments must be cleaned using one of the two techniques described below. It is recommended to start cleaning as soon as possible to prevent any soil drying on the instrument.

Only agents, solutions and detergents of proven efficacy must be used. The agents, solutions and detergents listed below are those used by the manufacturer during the validation of these instructions. The use of other products may lead to incomplete cleaning or premature wear of the instrument.

5.2.1 Manual cleaning

Equipment and consumables:

- Enzyme solution Alkazyne®
- Alkaline detergent Neodisher® MediClean Forte
- Osmosed water or water of equivalent chemical and microbiological quality
- Soft nylon bristle brush
- Disposable towels
- Ultrasonic bath

Method:

- 1) Rinse and brush the instrument for at least 30 seconds to remove excess soil while completely immersed in a bath of osmosed water or water of equivalent chemical and microbiological

quality at ambient temperature (15 to 25°C).

- 2) Completely immerse the instrument in the Alkazyne® 0.5% v/v enzyme solution (prepared according to the manufacturer's instructions) for 15 minutes at ambient temperature (15 to 25°C, ideally 20°C).
- 3) Brush all surfaces with a soft nylon bristle brush for at least 30 seconds until all visible soil is eliminated. Ensure that the grooves and holes are thoroughly cleaned.
- 4) Remove the instrument from the solution.
- 5) Rinse the instrument under running water at ambient temperature (15 to 25°C) for at least 3 minutes.
- 6) Perform ultrasonic cleaning of the instrument completely immersed in MediClean 0.5% v/v detergent (prepared according to the manufacturer's instructions) for 15 minutes at a maximum temperature of 40°C.
- 7) Rinse the instrument with osmosed water at ambient temperature (15 to 25°C) for at least 3 minutes. Check that the water enters the grooves and that the blind holes are filled and emptied several times.
- 8) Carefully dry the instrument with disposable towels.
- 9) Visually inspect the instrument in a well-lit room to confirm that there is no remaining soil.
- 10) Repeat steps 1 to 9 described above if any visible soil persists.

5.2.2 Automatic cleaning

Equipment and consumables:

- Enzyme solution Alkazyne®
- Alkaline detergent Neodisher® MediClean Forte
- Osmosed water and purified water or water of

- equivalent chemical and microbiological quality
- Soft nylon bristle brush
- Ultrasonic bath
- Washer-disinfector validated and maintained in compliance with current local procedures

Method:

- 1) Rinse and brush the instrument for at least 30 seconds to remove excess soil while completely immersed in a bath of osmosed water or water of equivalent chemical and microbiological quality at ambient temperature (15 to 25°C).
- 2) Perform ultrasonic cleaning of the instrument completely immersed in the Alkazyme® 0.5% v/v enzyme solution (prepared according to the manufacturer's instructions) for 15 minutes at ambient temperature (15 to 25°C).
- 3) Brush all surfaces with a soft nylon bristle brush for at least 30 seconds in the solution until all visible soil is eliminated. Ensure that the grooves and holes are thoroughly cleaned.
- 4) Rinse the instrument carefully with purified water at ambient temperature (15 to 25°C) for at least 1 minute.
- 5) Load the instrument into the washer-disinfector, placing it in a position to ensure drainage of the grooves and holes.
- 6) Run the washer-disinfector for one 10-minute cycle at a temperature of 55°C with MediClean 0.5% v/v detergent (prepared according to the manufacturer's instructions).
- 7) When unloading, visually inspect the grooves, the holes and other difficult to access zones of the instrument in a well-lit room to check that all visible soil has been eliminated. If necessary, repeat the cycle and/or clean manually.
- 8) In addition, check that the instrument is com-

pletely dry. If necessary, use disposable towels to remove any possible traces of water.

5.3 Inspection

Visually inspect the instruments in a well-lit room to detect any sign of corrosion, damage and wear. The cutting edges must be uniform and free from nicks.

Remove any damaged or blunt instruments, cleaning them of all biological substances and place them for disposal in compliance with current laws and regulations.

5.4 Sterilisation

5.4.1 Packaging

The instruments must be packed individually in a double medical quality sterilisation pouch. Ensure that the size of the packaging is sufficient so that the packed instrument does not strain the seams or tear the pouch. Also ensure as far as possible that the instruments do not collide with each other.

5.4.2 Cycles

The instruments must be moist heat sterilised in an autoclave with an active air removal cycle (ISO 17665-1). The autoclave must be validated, maintained and calibrated in compliance with current local procedures.

The following cycles have been validated to ensure a sterility assurance level (SAL) of 10^{-6} :

Type of cycle	Active air removal		
Temperature	132°C	134°C	134°C
Exposure time (minimum)	4 minutes	3 minutes	18 minutes*
Drying time (minimum)	30 minutes	30 minutes	30 minutes

* Moist air sterilisation parameters recommended by the World Health Organisation (WHO) for instruments for which there is a risk of TSE/CJD (transmissible spongiform encephalopathy and Creutzfeldt-Jakob disease) contamination.

6. STORAGE AND DISPOSAL

Instruments at the end of their life must be stored protected from dust, moulds, insects and pests avoiding extremes of temperatures and humidity.

When the inspection step (cf. chapter 5.3) designates an instrument as unsuitable for use, it must be cleaned of all biological substances and placed for disposal in compliance with current laws and regulations.

A shelf-life for sterilised instruments must be defined by each healthcare centre.

7. DEFINITION OF SYMBOLS

Symbol	Description
 2862	CE marking with identification number of the notified body.
	Manufacturer
	Consult the instructions for use
	Non sterile
	Authorised representative in the European Union
	Lot code
	Catalogue reference
	Unique device identifier
	Quantity
	Medical Device

8. SERIOUS INCIDENTS

In case of a serious incident related to the instrument, the user is responsible for reporting this event to the manufacturer, to the competent authority of the Member State in which the user is established, to the distributor and if appropriate, to the representative.

The contact details of the manufacturer and the European representative are available on the back cover.

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MPS Micro Precision Systems AG
Emergo Europe B.V.