MÈCHES – Notice d'utilisation	FR
СВРЕДЛАТА – Инструкции за употреба	BG
VRTACÍCH NÁSTAVCŮ – Návod k použití	CS
BOREKRONER – Brugsanvisning	DA
BOHRER – Gebrauchsanleitung	DE
ΤΡΥΠΑΝΙΩΝ – Οδηγίες Χρήσης	EL
DRILL BITS – Instructions for use	EN
BROCAS – Manual de uso	ES
PUURIOTSAKUD – Kasutusjuhend	ET
FÚRÓFEJEK – Használati utasítás	HU
PUNTE PER TRAPANO – Istruzioni per l'uso	IT
GRAŽTŲ – Naudojimo instrukcija	LT
URBJU ASMEŅUS – Lietošanas instrukcija	LV
BOREKRONENE – Bruksanvisning	NO
WIERTŁA – Instrukcja użytkowania	PL
BROCAS – Manual de Utilização	PT
VÂRFURILOR DE BURGHIIU – Instrucțiuni de utilizare	RO
VRTÁKOV – Návod na použitie	SK
SVEDROV – Navodila za uporabo	SL
BORRKRONORNA – Bruksanvisning	SV
MATKAP – Kullanım Talimatları	TR

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DRILL BITS – Instructions for use

1. INTRODUCTION	47
2. OVERVIEW	47
3. WARNINGS AND PRECAUTIONS FOR USE	48
4. INSTRUCTIONS FOR USE	49
5. REPROCESSING	49
6. STORAGE AND DISPOSAL	51
7. DEFINITION OF SYMBOLS	52
8. SERIOUS INCIDENTS	52

EN

1. INTRODUCTION

This document provides instructions for use and reprocessing for users of the drill bits marketed under the name of the manufacturer.

It contains the instructions for use of the instruments and also instructions for cleaning, inspection, sterilisation and storage.

The cleaning and sterilisation procedures described in this manual have been validated by the manufacturer. Other reprocessing methods may be suitable, but they must be validated beforehand by the end user and remain the responsibility of the latter. Moreover the end user must comply with the laws and regulations in countries with stricter reprocessing requirements than those specified in this manual.

2. OVERVIEW

2.1 Description and composition

Drill bits are reusable instruments designed to drill cylindrical holes in bone for various surgical procedures.

Drill bits are composed of stainless steel.

The main action is obtained by cutting edges, and their flutes allow debris to be removed.

2.2 Intended use

The instruments described in this manual are used to drill cylindrical holes in bone at a precise location through a guide via a tool holder or directly connected to a surgical motor providing rotational power.

The instruments have been validated for 20 cycles of use, including the required reprocessing (cleaning and sterilisation), without encountering any functional problem. Any further use is the respon-

sibility of the user. The use of the instrument is not limited by the number of validation cycles, but must only be reused if inspection after reprocessing does not show any signs of wear or damage such as described in [chapter 5.3](#).

They must only be used by qualified personnel fully trained in their use and in the appropriate surgical procedures in an operating theatre.

2.3 Contra-indications

Drill bits are not intended for use on children because reaming the growth plate impairs bone growth. The use of drill bits in pregnant women is also forbidden because of the contra-indication to anaesthesia in this patient population.

Drill bits are not intended for use in maintenance operations such as for sharpening cutting edges.

2.4 Performance characteristics

The cutting diameter of drill bits is from 1.5 to 4.5mm and a cutting depth of 15 to 50mm. They are designed for a rotation speed of 150 to 1200 RPM.

3. WARNINGS AND PRECAUTIONS FOR USE

3.1 Prerequisites before use

A full cycle of reprocessing consisting of cleaning, inspection and sterilisation must be performed before each use and this applies to new instruments as well. The use of non-reprocessed instruments can damage the patient's health by causing infection.

An inspection must be performed before each use. Any instrument showing signs of wear or damage must under no circumstances be used and must be

replaced by a new instrument.

Before use, the drill bits, the connections between the various parts, including that with the surgical motor, must be duly checked. If the coupling between the parts breaks during surgery, the duration of the procedure will be extended.

3.2 Precautions for use

The use of personal protective equipment is strongly recommended during the handling of contaminated or potentially contaminated instruments.

Careful handling of the drill bits is recommended to avoid the risk of injury and tearing of surgical gloves on the cutting edges.

In case of contact between instruments and any other equipment or tools present, specifically within the location during use, metal particles may be generated in the surgical cavity due to friction.

It is recommended to start the power tool gradually before applying axial force. This allows early detection and reaction in case of interference and reduce the risk of applying an excessive torque and/or an abrupt stoppage of the instrument. In case of resistance, do not force the instrument, but reverse the drive direction and remove it from the cavity.

When drilling obliquely to the surface of the bone, the point of the drill can easily slide along the bone. It is helpful first to orient the drill perpendicular to the near cortex to start the hole before gradually reorienting it to the needed angle.

3.3 Environment for use

The products covered by this manual are intended for use in operating theatres in sterile areas, in the human body including contact with vital fluids such as blood.

During the entire life cycle of the product, the maximum permissible temperature for instruments is 137°C, beyond which they may deteriorate. In addition, strongly alkaline solutions (pH > 11) and hypochlorite solutions should be avoided because they promote corrosion of metallic parts.

3.4 Reprocessing precautions

During cleaning processes, personal protective equipment usage is strongly recommended because contact with instruments soiled with human blood may cause infection and contamination of staff.

It is strictly prohibited to use metal brushes for cleaning as they cause premature wear of the instrument.

Drill bits are simple instruments that require no special attention during cleaning.

The cleaning process of the drill bits must start as soon as possible after use, soiled instruments that have dried are more difficult to clean.

As far as possible, any abrupt contact with the cutting parts of the drill bits must be avoided because there is a risk of injury and tearing of surgical gloves.

4. INSTRUCTIONS FOR USE

Before use, consult the precautions for use described in [chapter 3.2](#).

The good condition of the drill bits must be checked before use.

All instruments should be used and connected to a surgical motor with the appropriate connector for the type of drill bit, in a sterile area of the operating theatre.

The bone drilling procedure shall be performed in compliance with current surgical techniques.

Before reaming ensure that the bone is sufficiently accessible.

5. REPROCESSING

Before reprocessing, consult the reprocessing precautions described in [chapter 3.4](#).

5.1 Reprocessing at the point of use

It is recommended to remove excess soil with lint-free disposable wipes as promptly as possible after the end of the procedure.

Transport the instruments wrapped in damp paper between the various sites of use and/or reprocessing, avoiding shocks and mechanical damage.

5.2 Cleaning

The instruments must be cleaned using one of the two techniques described below. It is recommended to start cleaning as soon as possible to prevent any soil drying on the instrument.

Only agents, solutions and detergents of proven efficacy must be used. The agents, solutions and detergents listed below are those used by the manufacturer during the validation of these instructions. The use of other products may lead to incomplete cleaning or premature wear of the instrument.

5.2.1 Manual cleaning

Equipment and consumables:

- Enzyme solution Alkazyme®
- Alkaline detergent Neodisher® MediClean Forte
- Osmosed water or water of equivalent chemical

and microbiological quality

- Soft nylon bristle brush
- Disposable towels
- Ultrasonic bath

Method:

- 1) Rinse and brush the instrument for at least 30 seconds to remove excess soil while completely immersed in a bath of osmosed water or water of equivalent chemical and microbiological quality at ambient temperature (15 to 25°C).
- 2) Completely immerse the instrument in the Alkazyne® 0.5% v/v enzyme solution (prepared according to the manufacturer's instructions) for 15 minutes at ambient temperature (15 to 25°C, ideally 20°C).
- 3) Brush all surfaces with a soft nylon bristle brush for at least 30 seconds until all visible soil is eliminated. Ensure that the grooves and holes are thoroughly cleaned.
- 4) Remove the instrument from the solution.
- 5) Rinse the instrument under running water at ambient temperature (15 to 25°C) for at least 3 minutes.
- 6) Perform ultrasonic cleaning of the instrument completely immersed in MediClean 0.5% v/v detergent (prepared according to the manufacturer's instructions) for 15 minutes at a maximum temperature of 40°C.
- 7) Rinse the instrument with osmosed water at ambient temperature (15 to 25°C) for at least 3 minutes. Check that the water enters the grooves and that the blind holes are filled and emptied several times.
- 8) Carefully dry the instrument with disposable towels.
- 9) Visually inspect the instrument in a well-lit room to confirm that there is no remaining soil.
- 10) Repeat steps 1 to 9 described above if any visible soil persists.

5.2.2 Automatic cleaning

Equipment and consumables:

- Enzyme solution Alkazyne®
- Alkaline detergent Neodisher® MediClean Forte
- Osmosed water and purified water or water of equivalent chemical and microbiological quality
- Soft nylon bristle brush
- Ultrasonic bath
- Washer-disinfector validated and maintained in compliance with current local procedures

Method:

- 1) Rinse and brush the instrument for at least 30 seconds to remove excess soil while completely immersed in a bath of osmosed water or water of equivalent chemical and microbiological quality at ambient temperature (15 to 25°C).
- 2) Perform ultrasonic cleaning of the instrument completely immersed in the Alkazyne® 0.5% v/v enzyme solution (prepared according to the manufacturer's instructions) for 15 minutes at ambient temperature (15 to 25°C).
- 3) Brush all surfaces with a soft nylon bristle brush for at least 30 seconds in the solution until all visible soil is eliminated. Ensure that the grooves and holes are thoroughly cleaned.
- 4) Rinse the instrument carefully with purified water at ambient temperature (15 to 25°C) for at least 1 minute.
- 5) Load the instrument into the washer-disinfector, placing it in a position to ensure drainage of the

grooves and holes.

- 6) Run the washer-disinfector for one 10-minute cycle at a temperature of 55°C with MediClean 0.5% v/v detergent (prepared according to the manufacturer's instructions).
- 7) When unloading, visually inspect the grooves, the holes and other difficult to access zones of the instrument in a well-lit room to check that all visible soil has been eliminated. If necessary, repeat the cycle and/or clean manually.
- 8) In addition, check that the instrument is completely dry. If necessary, use disposable towels to remove any possible traces of water.

5.3 Inspection

Visually inspect the instruments in a well-lit room to detect any sign of corrosion, damage and wear. The cutting edges must be uniform and free from nicks.

Remove any damaged or blunt instruments, cleaning them of all biological substances and place them for disposal in compliance with current laws and regulations.

5.4 Sterilisation

5.4.1 Packaging

The instruments must be double-wrapped individually in a medical quality sterilisation pouch. Ensure that the size of the packaging is sufficient so that the wrapped instrument does not strain the seams or tear the pouch. Also ensure as far as possible that the instruments do not collide with each other.

5.4.2 Cycles

The instruments must be moist heat sterilised in an autoclave with an active air removal cycle (ISO 17665-1). The autoclave must be validated, main-

tained and calibrated in compliance with current local procedures.

The following cycles have been validated to ensure a sterility assurance level (SAL) of 10⁻⁶:

Type of cycle	Active air removal		
Temperature	132°C	134°C	134°C
Exposure time (minimum)	4 minutes	3 minutes	18 minutes*
Drying time (minimum)	30 minutes	30 minutes	30 minutes

* Moist air sterilisation parameters recommended by the World Health Organisation (WHO) for instruments for which there is a risk of TSE/CJD (transmissible spongiform encephalopathy and Creutzfeldt-Jakob disease) contamination.

6. STORAGE AND DISPOSAL

Instruments at the end of their life must be stored protected from dust, moulds, insects and pests avoiding extremes of temperatures and humidity.

When the inspection step (cf. [chapter 5.3](#)) designates an instrument as unsuitable for use, it must be cleaned of all biological substances and placed for disposal in compliance with current laws and regulations.

A shelf-life for sterilised instruments must be defined by each healthcare centre.

7. DEFINITION OF SYMBOLS

Symbol	Description
	CE marking with identification number of the notified body.
	Manufacturer
	Consult the instructions for use
	Non sterile
	Authorised representative in the European Union (mandataire)
	Lot code
	Catalogue reference
	Unique device identifier
	Quantity
	Medical Device

8. SERIOUS INCIDENTS

In case of a serious incident related to the instrument, the user is responsible for reporting this event to the manufacturer, to the competent authority of the Member State in which the user is established, to the distributor and if appropriate, to the representative.

The contact details of the manufacturer and the European representative are available on the back cover.

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MPS Micro Precision Systems AG
Emergo Europe B.V.