

Let's innovate together

m^ps  PRECIMED



Broach Handle


General documentation

[Technical Data](#)


www.mpsprecimed.com

Broach Handle


Customized to your requirements



Model A
Double offset Broach Handle




Model B
Double offset Broach Handle




Model C
Simple offset Broach Handle



Model D
Straight Broach Handle



Model E
Straight Broach Handle



Model F
Double offset Broach Handle

Model	Designation	Surgical approach	Offset 1	Offset 2	Length	Compatibility with interface MPS
A	Double offset Broach Handle (Must be ordered as a pair)	Antero-lateral (Roettinger)	Recommended : 12 mm	Fixed at 28mm (Preformed raw material)	380 mm max.	Male : Yes Female : Yes
B	Double offset Broach Handle (Must be ordered as a pair)	Anterior	Recommended Range : 30 to 70 mm	Fixed at 43.5mm (Preformed raw material)	280 mm max.	Male : Yes Female : Yes
C	Simple offset Broach Handle	Anterior	Recommended Range : 40 to 90 mm	-	Recommended : between 250 and 300 mm	Male : Yes Female : Yes
D	Straight offset Broach Handle	Posterior / Lateral	-	-	Recommended : between 250 and 300 mm	Male : Yes Female : Yes
E	Straight offset Broach Handle	Anterior (Hueter)	-	-	Recommended : between 250 and 300 mm	Male : Yes Female : Yes
F	Double offset Broach Handle (Must be ordered as a pair)	Anterior (Hueter)	Recommended : 12 mm	Fixed at 28mm (Preformed raw material)	380 mm max.	Male : Yes Female : Yes

In all cases MPS Precimed will adapt the interface to suit the customer requirements.

Note: Impaction plate can be orientated per request. Dimensions are fixed (50 x 80 mm).

